

DOC NO	61905
PAGE	1 of 6
REV	7.0

SPECIFICATION – PRINTED MATERIALS

Rev	Change Note	Date document drafted	Document prepared by (name)	Document Checked by (name)
3.0	C21846	21 Nov 08	Stuart Plascott	Sandra Lee
4.0	C22301	27 Jan 09	Andrew Hung	Sandra Lee
5.0	C23122	01 June 09	Andrew Hung	Sandra Lee
6.0	C23602	12 Aug 09	Andrew Hung	Sandra Lee
7.0	C24245	09 Nov 09	Andrew Hung	Sandra Lee

MASK FITTING TEMPLATES

**NOTE: MANUFACTURER TO PRINT FROM NATIVE FILES SUPPLIED BY RESMED
DO NOT PRINT FROM THIS PDF**

1. DETAILS

Style: Screenprinted film panel, die cut

OR

Laser cut

Tools: TBA

Overall Dimensions: **All cards have the same overall dimensions**

Height 175mm x **Width** 140mm

Tolerance: ± 0.5mm

Thickness: 0.8mm

Tolerance: ± 0.1mm

Refer to individual sheets for detailed dimensions.

Materials: 0.8 mm extruded sheet Polypropylene.

Opaque, glossy white, no texture

Food contact grade. Random co-polymer.

Alignment: Graphic misalignment not to exceed 0.25mm.

Artwork: **Page 1: 61905 Cover pg 1.eps**

Page 2: 61905 Quattro Liberty pg 2.eps

Page 3: 61905 Micro Act LT SoftGel pg 3.eps

Page 4: 61905 UM II Activa pg 4.eps

If retypeset, the same styles and sizes must be maintained.

Printing: All artwork should be printed on the top surface as shown below.


Colours: Additional colour is PMS 287

DOC NO	61905
PAGE	2 of 6
REV	7.0

SPECIFICATION – PRINTED MATERIALS

Assembly:

All four templates are to be assembled using a Loose Leaf Metal Ring similar to that shown below.




2. RESMED QC INSPECTION

Refer to AWI203-002 - Printed Materials/Labels/CDs/Videos Section 19
Agreed representative samples are to be held, one by ResMed and one by the manufacturer.
Each shall be identified by the ResMed "Representative Sample Label".


DOC NO	61905
PAGE	3 of 6
REV	7.0

SPECIFICATION – PRINTED MATERIALS




DOC NO	61905
PAGE	4 of 6
REV	7.0

SPECIFICATION – PRINTED MATERIALS



DOC NO	61905
PAGE	5 of 6
REV	7.0

SPECIFICATION – PRINTED MATERIALS



DOC NO	61905
PAGE	6 of 6
REV	7.0

SPECIFICATION – PRINTED MATERIALS

